

# Quatromyicide II

## Disinfectant & Sanitizer



### Applications:

Dilute according to the dilution rates mentioned on the label in potable water. Apply using a mop, sponge, cloth or trigger sprayer.

**Disinfection & Sanitizing:** This product is an effective sanitizer and disinfectant, when used at the dilution rates indicated, on visibly clean, hard, nonporous, inanimate surfaces. As this product's dilution rate and procedure differs according to its application, please refer to the label for the correct usage for this product.

### Benefits:

- Disinfects, sanitizes, and deodorizes
- Effectively removes mold and mildew
- Versatile – dilute according to your cleaning application
- Ideal for the food industry and delicate personnel settings such as daycares or senior care facilities
- Suitable for all washable hard surfaces
- Available as concentrate or ready-to-use

### Codes:

- **Product Code:** 52891
- **UPC:** 0 56069 52891 7

### Specifications:

Colour	Colourless
DIN	02321963
Contact Time	10 minutes
Odour	Mild
pH	7.5 - 8.0
Viscosity	Not determined
Solvents	No
Phosphates	No
WHMIS	Exempt
ECOLOGO®	No
Food Safe	Yes
Kosher	Yes
Shelf Life	Expiry Date

### Dilution Rates:


\*Tip colour proportioning systems at 1 GPM / 3.5 GPM

#### Third Sink Sanitizing

60 ml / 30 L

#### Disinfection

1:250



### Available Format:

\* Refer to Quatromyicide II RTU Spec Sheet



750 ml



5 L

### Packaging:

**Bottle:** 6.1" x 6.1" x 12.6" (15.5 x 15.5 x 31.9 cm)  
12.25 lbs (5.56 kg)

**Box:** 12.4" x 12.4" x 13"  
49 lbs (22.23 kg)  
4 per case

